


1 MALE ADAPTERS




Rotating male luer lock for tube \varnothing 6.8 mm
Material: Rigid PVC - PC / STD Pack: 10.000 pcs

- A00109990** - ETO (Neutral)
- A00109992** - ETO (Red)
- A00109993** - ETO (Blue)
- A00109980** - Radiation (Neutral)
- A00109982** - Radiation (Red)
- A00109983** - Radiation (Blue)




Rotating male luer lock for tube \varnothing 5.5 mm
Material: Rigid PVC - PC / STD Pack: 10.000 pcs

- A00119990** - ETO (Neutral)
- A00119992** - ETO (Red)
- A00119993** - ETO (Blue)
- A00119980** - Radiation (Neutral)
- A00119982** - Radiation (Red)
- A00119983** - Radiation (Blue)



Rotating male luer lock for tube \varnothing 4.1 mm
Material: Rigid PVC - PC / STD Pack: 10.000 pcs


- A00129990** - ETO (Neutral)
- A00129992** - ETO (Red)
- A00129993** - ETO (Blue)
- A00129980** - Radiation (Neutral)
- A00129982** - Radiation (Red)
- A00129983** - Radiation (Blue)



1.2 - Fixed male luer lock with vented cap


Male luer lock with vented cap for tube Ø 2.0 mm
Material: ABS - LDPE / STD Pack: 10.000 pcs

A01409990 - ETO - Radiation



Male luer lock with vented cap for tube Ø 2.5 mm
Material: ABS - LDPE / STD Pack: 10.000 pcs

A01419990 - ETO - Radiation



Male luer lock with vented cap for tube Ø 3.2 mm
Material: ABS - LDPE / STD Pack: 10.000 pcs

A01479990 - ETO - Radiation



Male luer lock with vented cap for tube Ø 3.7 mm
Material: ABS - LDPE / STD Pack: 10.000 pcs

A01489990 - ETO - Radiation



Male luer lock with vented cap for tube Ø 4.1 mm
Material: ABS - LDPE / STD Pack: 10.000 pcs

A01429990 - ETO - Radiation



Male luer lock with vented cap for tube Ø 5.5 mm
Material: ABS - LDPE / STD Pack: 10.000 pcs

A01449990 - ETO - Radiation




Male luer lock with vented cap for tube Ø 6.8 mm
Material: ABS - LDPE / STD Pack: 10.000 pcs

A01459990 - ETO - Radiation




Male luer lock for tube Ø 2.0 mm
Material: ABS / STD Pack: 20.000 pcs

C13304020 - ETO - Radiation




Male luer lock for tube Ø 2.5 mm
Material: ABS / STD Pack: 20.000 pcs

C13314020 - ETO - Radiation



Male luer lock for tube Ø 3.2 mm
Material: ABS / STD Pack: 20.000 pcs

C13374020 - ETO - Radiation



Male luer lock for tube Ø 3.7 mm
Material: ABS / STD Pack: 20.000 pcs

C13384020 - ETO - Radiation



Male luer lock for tube Ø 4.1 mm
Material: ABS / STD Pack: 20.000 pcs

C13324020 - ETO - Radiation



Male luer lock for tube Ø 5.5 mm
Material: ABS / STD Pack: 20.000 pcs

C13344020 - ETO - Radiation



Male luer lock for tube Ø 6.8 mm
Material: ABS / STD Pack: 20.000 pcs

C13354020 - ETO - Radiation

