


# 4

## CLAMPING DEVICES




Small clamp on/off for tubes Ø 4.1 mm / Ø 5.5 mm  
Material: PP / STD Pack: 8.000 pcs

**C13203131** - ETO - Radiation (Neutral)  
**C13203132** - ETO - Radiation (Red)  
**C13203133** - ETO - Radiation (Blue)




Medium clamp on/off for tubes Ø 5.5 mm / Ø 6.8 mm  
Material: PP / STD Pack: 3.000 pcs

**C13503131** - ETO - Radiation (Neutral)  
**C13503132** - ETO - Radiation (Red)  
**C13503133** - ETO - Radiation (Blue)




Big clamp with closure regulation for tube Ø 6.8 mm  
Material: PP / STD Pack: 2.000 pcs

**C20213030** - ETO - Radiation (Neutral)  
**C20213032** - ETO - Radiation (Red)  
**C20213033** - ETO - Radiation (Blue)




Big clamp on/off for tube Ø 6.8 mm  
Material: PP / STD Pack: 2.000 pcs

**C20303070** - ETO (Neutral)  
**C20303072** - ETO (Red)  
**C20303073** - ETO (Blue)




Adjustable hose clamp, 1.5 x 12 mm  
length 355 mm  
Material: LDPE / STD Pack: 4.000 pcs

**B90003B** - ETO - Radiation (Blue)



Roller housing for tube Ø 4.1 mm  
Material: ABS / STD Pack: 3.500 pcs

**C13004061** - ETO - Radiation (White)



Roller housing wheel for tube Ø 4.1 mm  
Material: ABS / STD Pack: 20.000 pcs

**C13104061** - ETO - Radiation (White)



Small roller housing for tube Ø 4.1 mm  
Material: PP / STD Pack: 12.500 pcs

**C28503021** - ETO (White)




Small roller housing wheel for tube Ø 4.1 mm  
Material: PP / STD Pack: 50.000 pcs

**C28603022** - ETO (Red)




Mini heparine clip  
STD Pack: 25.000 pcs

**C12501030** - ETO (Neutral) Rigid PVC  
**C12502120** - ETO - Radiation (Neutral) HDPE




Medium heparine clip  
STD Pack: 25.000 pcs

**C12601030** - ETO (Neutral) PVC  
**C12602120** - ETO - Radiation (Neutral) HDPE



Out-distance clip 2 channels for tube Ø 3.1 mm  
Material: PP / STD Pack: 10.000

**C18003030** - ETO - Radiation (Neutral)




Out-distance clip 2 channels for tube Ø 4.1 mm  
Material: PP / STD Pack: 10.000

**C18013030** - ETO - Radiation (Neutral)




Lam flow regulator "0 / 250 ml" tube Ø 4.1 mm  
for gravity infusion  
Material: ABS - Synthetic Rubber / STD Pack: 1.000 pcs

**A01509991** - ETO - Radiation (white)




Clamp sheet for bag  
Material: ABS / STD Pack: 4.000 pcs

**B10004B** - ETO - Radiation (Blue)



Open slide clamp for tubes  
Material: ABS / STD Pack: 5.000 pcs

**B10003B** - ETO - Radiation (Blue)



Slide clamp for tubes  
Material: ABS / STD Pack: 5.000 pcs

**B10029B** - ETO - Radiation (Blue)

